

# *Relieva* **Ultirra**<sup>™</sup>

Sinus Balloon Catheter

- > EASY PROCEDURAL IRRIGATION
- > ENHANCED BALLOON CATHETER PERFORMANCE
- > DURABLE SINUS DILATION

Acclarent<sup>®</sup>


# Relieva Ultirra™ Sinus Balloon Catheter

The *Relieva Ultirra* Sinus Balloon Catheter combines high performance balloon dilation with intra-sinus irrigation to deliver clinically significant improvement in sinus symptoms



## EASY PROCEDURAL IRRIGATION IS CONVENIENT IN EVERY CASE

- > Displace sinus contents using a combination of side and distal-tip jets
- > Increase procedural efficiency by eliminating device exchanges
- > Target irrigation by positioning the soft catheter tip deep within the sinus



## ENHANCED BALLOON CATHETER PROVIDES SUPERIOR PERFORMANCE

- > Smoother passage through sinus anatomy\*
- > Easier insertion and retraction through the sinus guide catheter\*
- > Success through multiple sinus dilations with robust balloon catheter performance



## DURABLE SINUS DILATION IN SIX BALLOON SIZES

- > Select dilation diameter and length to address specific sinus anatomy needs
- > Clinically proven to restore sinus drainage and function with minimal tissue disruption
- > Safe and effective dilation in pediatric maxillary sinuses

Catalog Number	Balloon Diameter (mm)	Balloon Length (mm)	Units (per box)
BC3512RU	3.5	12	5
BC0516RU	5.0	16	5
BC0616RU	6.0	16	5
BC0716RU	7.0	16	5
BC0524RU	5.0	24	5
BC0724RU	7.0	24	5

\*As compared to *Relieva Solo Pro*™. Data on file.

For more information on *Relieva Ultirra* and Acclarent, Inc., please contact your Acclarent representative or visit [www.acclarent.com](http://www.acclarent.com).

*Balloon Sinuplasty*™ System is intended for use by or under the direction of a physician that is trained in the use of *Balloon Sinuplasty* Technology. Prior to use, it is important to read the Instructions for Use and to understand the contraindications, warnings, and precautions associated with this device.

1525-B O'Brien Drive  
Menlo Park, CA 94025  
U.S.A.

t 650.687.5888  
f 650.687.5886  
1.877.775.2789



[www.acclarent.com](http://www.acclarent.com)

[info@acclarent.com](mailto:info@acclarent.com)

© 2011 Acclarent, Inc. All Right Reserved. MKT02194 Rev. B